

NEMA PET Sensitivity Phantom (NU 2-2018)[™]

Model PET/NEMA-SEN/P

Main Features

- 6 Concentric aluminum tubes used to detect camera sensitivity in PET

Specifications

5 internally stacked aluminum tubes all 700 mm in length

1st Tube inside diameter 3.9 mm, outside diameter 6.4 mm

2nd Tube inside diameter 7.0 mm, outside diameter 9.5 mm

3rd Tube inside diameter 10.2 mm, outside diameter 12.7 mm

4th Tube inside diameter 13.4 mm, outside diameter 15.9 mm

5th Tube inside diameter 16.6 mm, outside diameter 19.1 mm

The innermost tube, a fillable polyethylene tubing has an inside diameter of 1 mm, outside diameter 3 mm


Shipping

Carton: 6" x 6" x 36"

Weight: 3lbs.



Close up of end of NEMA PET Sensitivity Phantom (NU 2-2018)[™]



Set of aluminum tubes used in NEMA PET Sensitivity Phantom (NU 2-2018)[™]



Tubing kit for NEMA Sensitivity PET Phantom (NU 2-2018)[™]



NEMA Sensitivity PET Phantom (NU 2-2018)[™]