


Solid Sphere Set (6)™



Solid Sphere Set (6)™

Solid Sphere Set™

Model ECT/SS/SET6

Main Features:

- Designed for use in all Data Spectrum circular and elliptical ECT phantoms

Main Applications:

- Simulates cold spherical “lesions”
- Quantitative evaluation of spatial resolution/object size, attenuation and scatter effects
- Evaluation of quantitative ECT reconstruction methods
- Research

Specifications:

Diameters of solid spheres: 9.5, 12.7, 15.9, 19.1, 25. and 31.8 mm