



Dynamic Cardiac Phantom

Call for Price* – For North America Sales ONLY

Includes:

- a.) Custom torso phantom with liver insert, lung inserts and beating left ventricular insert with separate latex ventricular chamber and latex myocardial wall chamber. Each insert and chamber can be filled independently.
- b.) Three replacement sets of latex membranes.
- c.) Ventricular chamber pump, motor, cables and computer controlled software. Stroke volume and heart rate can be varied.
- d.) Ninety day warranty including parts and labor. User responsible for cost of shipping phantom to DSC should warranty service be required.
- e.) Assembly and User's manual

Excludes:

- a.) Required laptop (preferred) computer is not included with phantom. Laptop computer to be supplied by user. Computer must have at least one RS-232 serial port and windows XP™ operating system. A USB-to-RS232 converter will not function properly and cannot be used.

Caveats:

- a.) Software supplied as is. Software crashes may randomly occur. Software crashes are not covered by warranty.
- b.) Although motor-cooling fans are supplied, motor may overheat and protective circuitry may shut-down motor during certain rare operational conditions such as, for example, long acquisition times at large stroke volumes and/or high heart rates. After a period of cooling, motor will resume normal operation. Motor shut-down due to overheating is not covered by warranty.
- c.) Volume and shape of ventricle may possibly change during a run due to leaking between chambers. Volume and shape changes of ventricle is not covered by warranty.

Quote is in US dollars.

Shipping charges are based on actual UPS shipping charges, if using other shipper shipping payment will be collect.

*All quotes good for 30 days, please note pricing can change.