

Fillable Spine Insert™

Model ECT/FIL-SPINE/I (user filled)

Spine Insert, With Liquid Bone™

Model ECT/BONE-SPINE/I (filled with liquid bone)

Main Features

- Designed to be used with the Elliptical Lung-Spine Insert™ (Model ECT/LUNG-SPINE/I) or Elliptical Lung-Spine Phantom™ (Model ECT/LUNG-SPINE/P) or Anthropomorphic Torso Phantom™ (Model ECT/TOR/P)

Main Applications

- Improved spine attenuation characteristic over Teflon® rod

Specifications

All clear material: PMMA
Outside height: ~ 190 mm
Inside height: ~ 152 mm
Outside diameter: ~ 45 mm
Inside diameter: ~ 38 mm
Volume: ~ 170 cc

Shipping

Carton: 7" x 7" x 3"
Weight: 2lbs.



Fillable Spine Insert™